

# Data Sheet



## Excel Life™ Medical Range, Common Plate Horizontal Single Powerpoint 15A – Antimicrobial\*

with Power Available LED and marked “RCD PROTECTED”

CAT NO. **EC787/15NCAWE, EC787/15NCABU, EC787/15NCARD**



Cat No	EC787/15NCAWE, EC787/15NCABU, EC787/15NCARD
Description	Excel Life™ Medical Range Common Plate – Horizontal Single Powerpoint 15A – Antimicrobial with Power Available LED and marked “RCD Protected”
Supply voltage	230 – 240V a.c.
Frequency	50Hz
Current	15A
Protection against solid bodies	IK04
IP rating	IP30
Power Available LED	Amber
Colour	White(WE), Blue(BU) or Red(RD)
Plate size	74mmW x 118mmH x 10mmD
Mounting centres	84mm
Operating temperature	5° to 35°C
Glass/tile cut out	56mm x 96mm
Standard	AS/NZS 3112

# Excel Life™ Medical Range, Common Plate Horizontal Single Powerpoint 15A – Antimicrobial\*

with Power Available LED and marked “RCD PROTECTED”



CAT NO. **EC787/15NCAWE, EC787/15NCABU, EC787/15NCARD**



- Chemical Resistant and Antimicrobial
- Internal safety shutters for extra safety
- Power available LED (Amber)
- 15A rated
- More fixing options – 4 corners
- Retained terminal screws
- Marked with “RCD protected”
- Colour White, Blue (UPS), Red (Essential)

*Antimicrobial Specifications			
<b>Maintenance</b>	Surface cleaning with cloth: Hospital-grade water-based cleaners only. Do not use any solvent based cleaner such as acetone, tar-remover or trichloroethylene.		
<b>Antimicrobial Properties</b>	Challenge Organism	Polycarbonate alloy in facia parts	Polycarbonate alloy in switch dolly
		Reduction Coefficient/Kill Rate	
<b>Note:</b> 5 Log = 99.999% 4 Log = 99.99% 3 Log = 99.9% 2 Log = 99%	<b>Escherichia Coli</b>	2.55 Log / 99.72%	6.72 Log / 99.9999%
	<b>Staphylococcus aureus</b>	2.70 Log / 99.79%	5.22 Log / 99.999%
	<p>Tests performed according to JIS Z 2801:2010 (Antimicrobial products – Test for antimicrobial activity and efficacy). Test method is equivalent to procedure detailed in ISO 22196:2011 (Measurement of antibacterial activity on plastics and other non-porous surfaces). Test results available upon request.</p> <p>Not to be used as an alternative or substitute for standard hygienic practices, procedures, or protocols.</p> <p>This device has a long-term antimicrobial effect. Lifespan is dependent upon cleaning procedures and the use of recommended cleaning products.</p> <p>Product manufactured with antimicrobial material based on Ag+ ions. This technology limits the proliferation of bacteria on our devices without increasing antibiotic resistance.</p>		

Single product		
Cat No	Dimensions	Barcode
EC787/15NCAWE	85mmW x 117mmH x 50mmD	9321001428190
EC787/15NCABU	85mmW x 117mmH x 50mmD	9321001428206
EC787/15NCARD	85mmW x 117mmH x 50mmD	9321001428213

Inner box (QTY 10)		
Cat No	Dimensions	Barcode
EC787/15NCAWE	175mmW x 165mmH x 125mmD	9321001429005
EC787/15NCABU	175mmW x 165mmH x 125mmD	9321001429012
EC787/15NCARD	175mmW x 165mmH x 125mmD	9321001429029

Outer box (QTY 120)		
Cat No	Dimensions	Barcode
EC787/15NCAWE	608mmW x 308mmH x 340mmD	9321001429814
EC787/15NCABU	608mmW x 308mmH x 340mmD	9321001429821
EC787/15NCARD	608mmW x 308mmH x 340mmD	9321001429838

HPM Legrand - AUS  
1300 369 777  
www.legrand.com.au  
ABN 31 000 102 661

HPM Legrand - NZ  
0800 476 009  
www.legrand.co.nz

Refer to your group buying office, HPM Legrand Sales Representative or HPM Legrand Sales Office for pricing.

Please note: specifications are subject to change without prior notice.